

## CRT-SM60G Glass Door Ambient Cabinet 60 litres

£895+VAT

Product Code: CRT-SM60G

### CRT-SM60G

Designed for the safe and correct storage of pharmaceuticals under limit regardless of surrounding room temperatures. This compact size space is a premium.

External Dimensions: 660H x 470W x 455D mm

Internal Dimensions: 525H x 375W x 310D mm

Capacity: Gross 60 litres / Net 57 litres

Suitable for Pharmaceutical Storage under +25°C MAX

Forced-Air Cooling for even temperature distribution

Digital LED Temperature Display with ShoreTouch controls

Audio / Visual Temperature Deviation Alarm

Automatic Defrosting

Door Lock and 2 keys

CFC/HFC & Ammonia free (R600a / 0GWP)

RPS, NHS compliant

Max. Ambient +32°C (Class SN)

220-240v 50/60 Hz

220-240v 50/60 Hz

5 Years Part & Labour Warranty (UK only)

5 Years Annual Calibration (UK only) on product registration

Certificate of Calibration included

