

## CRT-SM45 Solid Door Ambient Cabinet 45 litres

£695+VAT

Product Code: CRT-SM45

### CRT-SM45

Designed for the safe and correct storage of pharmaceuticals under limit regardless of surrounding room temperatures. This compact size space is a premium.

External Dimensions: 715H x 380W x 420D mm

Internal Dimensions: 445H x 300W x 250D mm

Capacity: Gross 45 litres / Net 42 litres

Suitable for Pharmaceutical Storage under +25°C MAX

Forced-Air Cooling for even temperature distribution

Digital LED Temperature Display with ShoreTouch controls

Audio / Visual Temperature Deviation Alarm

Automatic Defrosting

Door Lock and 2 keys

Adjustable Shelving x 2

CFC/HFC & Ammonia free (R600a / 0 GWP)

RPS, NHS, WHO, RCVS compliant

Max. Ambient +32°C (Class SN)

Made in England

5 Years Part & Labour Warranty (UK only)

5 Years Annual Calibration (UK only) on product registration

Certificate of Calibration included

