

SMP65 Portable Medical Refrigerator – 65 litre

£999+VAT

Product Code: SMP65 Portable

SMP65 -

Larger portable medical refrigerator (compressor type).

Designed to maintain +2°C to +8°C.

Ideal for mobile nurses and vets.

12 Volt DC or 100-240 Volt AC

External 450H x 725W x 510D mm

Capacity 65 litres

Temperature Range +2°C to +8°C

Forced-Air Cooling

Digital LED Temperature Display

Audio / Visual Temperature Deviation Alarm

MAX / MIN Temperature Memory

Removable Wire Basket

24 Kgs (two-man carry)

42 Watts

CFC & Ammonia free

Finished in durable-moulded exterior

Carry Handles

Max. Ambient +32°C

Assembled in England

2 Years Parts & Labour Warranty (UK only)

Supplied with Certificate of Calibration

