

SM1502 Solid Door Pharmacy Refrigerator – 151 litres

£1050+VAT

Product Code: SM1502

Ideal size for larger pharmacies or surgeries.

External Dimensions: 850H x 550W x 610D mm

Internal Dimensions: 720H x 455W x 440D mm (gross)

Capacity: Gross 151 litres / Net 148 litres

Temperature Range +2°C to +8°C

Suitable for Vaccine & Pharmaceutical Storage

Forced-Air Cooling for even temperature distribution

Digital LED Temperature Display with ShoreTouch controls

Audio / Visual Temperature Deviation Alarm

Integral MAX / MIN Temperature Memory

Automatic Defrosting

Door Lock & 2 keys

CFC & Ammonia free (R600a)

Adjustable Shelving x 3

RPS, NHS, WHO, RCVS compliant

220-230v 50Hz

Max. Ambient +32°C (Class SN)

Made in England

5 Years Part & Labour Warranty (UK only)

Certificate of Calibration included

