

SM161 Solid Door Pharmacy Refrigerator – 141 litres

£1150+VAT

Product Code: SM161

SM161 -

Compact under-counter size ideal for pharmacies and surgeries

External Dimensions: 835H x 600W x 615D mm

Internal Dimensions: 667H x 441W x 439D mm (gross)

Capacity: Gross 141 litres / Net 130 litres

Temperature Range +2°C to +8°C

Suitable for Vaccine & Pharmaceutical Storage

Forced-Air Cooling for even temperature distribution

Digital LED Temperature Display with ShoreTouch controls

Audio / Visual Temperature Deviation Alarm

Audio / Visual Temperature Deviation Alarm

Automatic Defrosting

Door Lock and 2 keys

High Efficiency - Thick Wall Insulation

Adjustable Shelving x 3

CFC & Ammonia free (R600a)

RPS, NHS, WHO, RCVS compliant

220-230v 50Hz

Max. Ambient +38°C (Class ST)

Made in England

5 Years Part & Labour Warranty (UK only)

Certificate of Calibration included

