

## SM45 Solid Door Pharmacy Refrigerator &#8211; 45 litre

£699+VAT

Product Code: SM45

Compact size ideal for nursing homes, and clinics.

External Dimensions: 715H x 380W x 420D mm

Internal Dimensions: 445H x 300W x 250D mm

Capacity: Gross 45 litres / Net 42 litres

Temperature Range +2°C to +8°C

Suitable for Vaccine & Pharmaceutical Storage

Forced-Air Cooling for even temperature distribution

Digital LED Temperature Display with ShoreTouch controls

Audio / Visual Temperature Deviation Alarm

Integral MAX / MIN Temperature Memory

Automatic Defrosting

Door Lock and 2 keys

CFC & Ammonia free (R600a)

Adjustable Shelving x 2

RPS, NHS, WHO, RCVS compliant

220-230v 50/60 Hz

Max. Ambient +32°C (Class SN)

Made in England

5 Years Part & Labour Warranty (UK only)

Certificate of Calibration included

